

HOME & INSTITUTIONAL MATTRESS GUIDE




All Specialty Product Orders or Inquiries Call
800-372-7749

www.blackburnsmed.com




THERAPEUTIC FOAM MATTRESSES



Therapeutic Foam Mattress

- 400# capacity
- Dual-layer blends patient comfort and clinical effectiveness
- Pressure redistribution across surface protects vulnerable skin
- Fluid-proof, wipe-clean bacteriostatically-treated top cover with hidden zipper



Therapeutic Foam Mattress

- 400# capacity
- Zoned construction with single layer value
- Bacteriostatically-treated, fluid-proof, wipe-clean cover with ultra-low moisture vapor transmission rate (MVTR)
- Hidden zipper design for easier cleaning and superior infection control



Therapeutic Raised Perimeter Mattress

- 400# capacity
- Raised perimeter features tissue-friendly Geo-Mattress® surface from edge-to-edge
- Helps prevent impingement of extremities
- Extra-firm edge bolster
- Shearing protection at transfer opening



Therapeutic Heavy Duty Solution for Bariatric Care Mattress

- 750# capacity
- The best of the clinically proven Geo-Mattress® series
- A heavy duty solution for bariatric care
- Fluid-proof bacteriostatically-treated, easily cleaned coated nylon cover

LOW-AIR-LOSS MATTRESSES

ALTERNATING LOW-AIR-LOSS



FUSION 2K™
Featuring SelectAir® Therapy

Preventative / Therapeutic Mattress

- 300# capacity
- Zoned, "True" low air loss system used to manage microclimate of the skin and improve moisture redistribution
- Length 80" or 84"; Width 36"
Fits standard hospital beds

SINGLE-THERAPY LOW-AIR-LOSS MATTRESS



SenTech
Stage IV 2000

- 500# weight capacity
- Single-zoned system
- Automatic pressures for weight and height settings
- Sensor technology to detect changes in patient's position and automatically adjust the airflow
- Mattress Widths 39", 42", 48", 54" & 60" available
- Mattress overlays available in 4" and 5" heights

BARIATRIC MULTI-THERAPY LOW-AIR-LOSS



SenTech
Stage IV 3000

- 500# capacity
- Multi-Zoned System offers three individual pressure zones for maximum patient comfort and therapy. Head, trunk, and foot sections can each be customized
- Automatic pressures for weight and height settings
- Sensor technology to detect changes in patient's position and automatically adjust the airflow
- Mattress Widths 39", 42", 48", 54" & 60" available